

# **FERNO Compact 1 & 2** Carry Chairs

Ideal for use in confined areas



Compact 1 Carry Chair



Compact 2 Carry Chair

**FOR MORE DETAILED  
INFORMATION SEE OVER**

## Specifications

	Height (folded)	Width (folded)	Depth (folded)	Weigh	Load Limit
Compact 1	935mm (765mm)	435mm (435mm)	430 mm (235mm)	9.1 kg	200kg/31st/440lbs
Compact 2	935mm (775mm)	445mm (445mm)	420 mm (285mm)	9.4 kg	200kg/31st/440lbs

**FERNO – is a world leading innovator and manufacturer of patient transport solutions.**

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Carry Chairs** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

**call** 01 5079900 or 069 66006  
**email** [sales@hayesfirstaid.com](mailto:sales@hayesfirstaid.com)  
**www** [www.hayesfirstaid.com](http://www.hayesfirstaid.com)



# FERNO Compact 1 & 2 Carry Chairs

The **Compact Carry Chair** range of chairs are durable, robust and are manufactured from stainless steel. Their one piece back and seat flame retardant, vinyl coated nylon covers allows the patient to adopt a comfortable position during transfer.

The **Compact range** design incorporates safety rings which slide over the folding mechanism.

- Seat covers can simply be wiped clean for maximum hygiene.
- Folds neatly into a compact size for easy stowage.
- 7" wheels for ease of manoeuvring.
- Maximum Load 200kg.
- Easily adjustable chest and ankle restraints are quickly applied
- Unique style lifting loops help the user to adopt a secure grip
- Large wheels and an axle bar to provide stability over uneven ground

The **Compact 1** proven chair is invaluable when transferring patients through corridors and lifts due to it's narrow design and large rear wheels.

The **Compact 2 Carry Chair** has been designed with ergonomic front grip handles for improved posture when carrying patients down stairs.



Compact 1 Carry Chair



Compact 2 Carry Chair



## Part Nos:

071481302/S	Compact 1 Carry Chair (black)
071481401/S	Compact 2 Carry Chair (black)
071481401/B10	Compact 2 Carry Chair (black) with Biosafe straps

## Accessories

038299401	Dust cover for Compact 1, & 2 (black)
081931400	Wall Bracket for Compact 1 & 2 chairs
031350600	Carry Case For Compact 1/2 Chair

FERNO Compact 1 & 2 Carry Chairs FPC002 10.14

**call** 01 5079900 or 069 66006  
**email** sales@hayesfirstaid.com  
**www** www.hayesfirstaid.com

Hayes First Aid, Castle Demesne, Newcastle West, Limerick, V42 N220, Ireland.

