

EVACUATION PRODUCTS

Evacuation products enable the safe and comfortable evacuation of mobility impaired people in the event of an emergency.

In the Evacusafe series there are several transit chairs, evacuation chairs and mattresses available.

EVACUSAFE Excel Evacuation Chair

The Excel Evacuation Chair is developed to evacuate mobility impaired people during the event of emergency.

The tracks create friction which acts as a brake to slow the descent and allows for smooth controlled operation to make sure the client will be safely transported to the place of safety.

This evacuation chair is suitable for a single person operation, where other evacuation products often require more help.



4 Wheel Transit Chair



Tri-Wheel Transit Chair



Evacuslider



EVACUSAFE Excel Evacuation Chair

The Excel Evacuation Chair is lightweight yet incredibly strong, the chair can transport persons of up to 182 kg, which makes it possible to transport children as well as heavy weighted people downstairs. With a supportive flat seat, armrests and padded headrest, the chair is very comfortable to stay in for extended periods.

Weight	13,6 kg
Maximum carrying weight	182 kg
Length folded	113 cm
Length operational	143 cm
Depth folded	29 cm
Depth operational	99 cm
Width	65 cm
To be operated by	1 person
Extra information	Padded headrest Flat, comfortable seat



EVACUSAFE 4 Wheel Transit Chair

The 4-Wheel Transit Chair is primarily used to transport persons across same level surfaces to a place of safety. This chair is similar to those used by hospitals and ambulance services. Basement areas can be tricky environments to evacuate from. In this situation the Tri-Wheel Transit Chair may be more suited to aid the evacuation.

Weight	9,6 kg
Maximum carrying weight	120 kg
Length folded	83 cm
Length operational	93 cm
Depth folded	26 cm
Depth operational	98 cm
Width	52 cm
To be operated by	2 to 4 persons
Extra information	Padded seat



EVACUSAFE Tri-Wheel Transit Chair

The Tri-wheel Transit Chair is developed to evacuate mobility impaired people in staircases. The unique 3 wheel system makes it easy to move on the steps and helps to evacuate the patient upstairs. The chair should be operated by 2 to 4 people, depending on the person who is being evacuated.

Weight	14 kg
Maximum carrying weight	90 kg
Length folded	91 cm
Length operational	88 cm
Depth folded	31 cm
Depth operational	128 cm
Width	53 cm
To be operated by	2 to 4 persons
Extra information	Flat, comfortable seat



EVACUSAFE Premium Evacuslider

The Premium Evacuslider is a cost efficient solution to safely evacuate mobility impaired people. The mattress is totally adjustable, which makes it possible to evacuate children as well as adults. It is made from a tough, wipe clean material and it provides comfort and support due to the thick foam mattress. The integrated foot pocket and 3 safety harnesses help to secure and cocoon the client. The flat backside makes it also possible to transport the patient horizontally.

Weight	5,5 kg
Maximum carrying weight	120 kg
Length folded	65 cm
Length operational	195 cm
Depth folded	25 cm
Depth operational	5 cm
Width	65 cm
To be operated by	2 persons
Extra information	Adjustable foot pocket