



Why choose ARKY from Hayes First Aid?

- ✓ Attractive price-quality ratio
- ✓ High quality materials and components
- ✓ All products are in stock, shipment within 24 hrs
- ✓ Cabinets fit defibrillators of all brands
- ✓ Customization/branding options
- ✓ All products are individually packed in protective plastic and a three-layer carton box

ARKY cabinets are the perfect solution for the safe and professional storage of an Automated External Defibrillator

METAL INDOOR CABINET

High quality

- Strong metal cabinet creates a seamless look
- Ideal for office buildings and public spaces
- Perfectly powder coated
- Logo engraved aluminum hand grip

Dimensions

- 38 x 38 x 20 cm
- Fits defibrillators of all brands

Options

- Acoustic alarm
- Available in green and white



Product no.	Description	Alarm	Fits defibrillators of all brands
60.111	White metal cabinet		✓
60.112	White metal cabinet with alarm	✓	✓
60.121	Green metal cabinet		✓
60.122	Green metal cabinet with alarm	✓	✓