

DefibStore 4000 AED Cabinet

Saving a life deserves the best protection

Overview

Hayes First Aid recognises the importance of a high-quality, secure cabinet to store your public access defibrillator. Manufactured from polycarbonate, the innovative DefibStore 4000 is the only polycarbonate AED cabinet on the market that can withstand high impact, is fire retardant to UL94 V0, and is UV stabilised against the effect of direct sun. A revolutionary cabinet that is uncompromising on quality and uncompromising on safety.



Key features

- Highly impact-resistant polycarbonate outdoor cabinet
- Fire retardant to UL94 V0 — will not support fire
- UV stabilised; will not be affected by sunlight
- Double-skin clamshell design with full IP66 rating
- Marine-grade stainless steel keypad lock
- Thermostatically controlled heater
- Battery-powered motion-sensor LED light, detachable as a torch
- Lightweight construction — only 10kg
- Fully recyclable material; 10-year limited guarantee

IP66

DUST & WATER TIGHT

UL94 V0

FIRE RETARDANT

10kg

LIGHTWEIGHT BUILD

10 YEAR

LIMITED GUARANTEE

Why the 4000 stands out

- **Unique double-skin construction** provides a robust, impact-resistant structure and creates additional insulation around the defibrillator.
- **Fully sealed clamshell design** — IP66 (EN60529) rated for a completely dust and water-tight enclosure, with a sloping top so surface water drains away from the door.
- **Marine-grade stainless steel keypad lock** salt-spray tested for 10,000 hours. Also available unlocked with a plastic winged handle and optional alarm.
- **Internal AED hook** elevates any defibrillator from the base and keeps it clearly visible through the large viewing window.
- **Simple to install** with a dual exit point for the power cable via the rear or base — no onsite drilling. Delivered fully assembled with plug, fixing kit and guide.

Order or enquire today

Email sales@hayesfirstaid.com Web www.hayesfirstaid.com

Made in the UK

MacroMed Healthcare Limited, trading as Hayes First Aid - Castle Demesne, Newcastle West, Co. Limerick, V42 N220, Ireland

Specifications

External dimensions	465mm (H) x 435mm (W) x 240mm (D)
Weight	10kg
Material	UV-stabilised polycarbonate
Ingress protection	IP66 (EN60529)
Fire rating	UL94 V0 (self-extinguishing)
Colours	High-vis yellow or green as standard
Heating & lighting	Thermostatic heater + motion-sensor LED light



How to order

Product code	Colour	Door locking	Heating & lighting	IP rating	Window
DFS4K442LLE-Y	Yellow	SS keypad lock	Yes	IP66	Yes
DFS4K442LLE-G	Green	SS keypad lock	Yes	IP66	Yes
DFS4K442ULE-Y	Yellow	Plastic wing handle	Yes	IP66	Yes
DFS4K442ULE-G	Green	Plastic wing handle	Yes	IP66	Yes

All cabinets are delivered fully assembled with plug, fixing kit, and installation & maintenance guide. The DefibStore 4000 is also available unlocked with optional alarm, and with a cooling facility — ask our sales team for details. Other colours available subject to minimum order quantities.

Order or enquire today

Email sales@hayesfirstaid.com Web www.hayesfirstaid.com

Made in the UK

MacroMed Healthcare Limited, trading as Hayes First Aid - Castle Demesne, Newcastle West, Co. Limerick, V42 N220, Ireland