

Sharpak 35 PLUS, 60 PLUS, 120 PLUS and 220 PLUS sharps containers

The Sharpak range of blue-top sharps containers are recommended for the segregation, collection and disposal of partially discharged sharps waste which are NOT contaminated with cytotoxic or cytostatic medicines.

This comprehensive range of 10 containers is available from 0.5L through to 26.3L and is accompanied by a range of brackets, trays and stands.

Sharpak sharps containers are approved to ISO 23907:2019 and tested and certified in accordance with the provisions of the United Nations recommendations of the United Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

Sharpak sharps containers MUST only be disposed by incineration.

Please note:- For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to the data sheet for the Code PURPLE range of Sharpak HYDRI sharps containers.

- MANUFACTURED, TESTED AND CERTIFIED TO ISO 23907:2019
- UN TYPE TESTED
- BLUE COLOUR CODED CLOSURES TO MEET THE RECOMMENDATIONS ON WASTE SEGREGATION AS DETAILED IN THE DEPARTMENT OF HEALTH'S TECHNICAL MEMORANDUM HTM-07-01
- SECURE CLOSURE WITH LOCKING FACILITY
- VISIOLOCK INDICATOR IS FITTED ON THE SHARPAK 35 PLUS, 60 PLUS, 120 PLUS, AND 220 PLUS,
- DISTINCTIVE COLOURING AND CLEAR LABELLING
- SEQUENTIAL NUMBERING FOR TRACEABILITY
- TOUGH AND DURABLE
- QUICK AND EASY TO ASSEMBLE
- COST EFFECTIVE
- PURPOSE DESIGNED KICK TROLLEY, WALLSTROLLEY BRACKETS, AND WALL UNIT
- INSULIN PEN NEEDLE REMOVING FEATURE ON SHARPAK 35 PLUS



PRODUCT	MAX CAPACITY	CAPACITY TO FILL LINE	PRODUCT CODE	QTY	HEIGHT	WIDTH	DEPTH
Sharpak 35 PLUS	3.5 litres	2.63 litres	SHK/CA1556HPB	40	195 mm	177 mm	206 mm
Sharpak 60 PLUS	6 litres	4.5 litres	SHK/CA1526HPB	25	229 mm	220 mm	234 mm
Sharpak 120 PLUS	12 litres	9.6 litres	SHK/CA1536HPB	15	325 mm	258 mm	268 mm
Sharpak 220 PLUS	22 litres	17.6 litres	SHK/CA1548HPB	8	383 mm	315 mm	327 mm
ACCESSORIES	CODE	QTY	ACCESSORIES	CODE	QTY		
Kick Trolley for 120 or 220 PLUS	SHK/SM16046	1	Wall & trolley bracket for 35 - 220 PLUS	SHK/SM05001	1		
Shaped procedure tray for 35 PLUS	SHK/SU03501	5	SAFE & SOUND wall unit for 60 PLUS	SHK/SS00060	1		

