

## TVAC

CE REF TM-909


NSN: 6515-01-540-7206

Vol. 20 per case, Total Weight: 4.8kg/10.5lbs

TyTek Medical's disposable portable suction unit with its resilient silicone bulb performs reliably in even the harshest climates without loss of suction. The TVAC Manual Suction Unit comes in a field friendly compact package, just pull out the built-in tube and start sucking. No bulky cartridges or brittle parts to break. After use, simply place the entire TVAC Manual Suction Unit back inside its packaging, seal, expose the biohazard label on the package and dispose appropriately. This entirely removes the possibility of cross contamination.

### Details:

- Strong, silicone bulb creates vacuum pressure of approx. 110 mmHg
- Collection bag volume >1,000ml
- Resilient, reliable silicone bulb prevents loss of suction in even the harshest climates
- Contoured suction tube tip prevents tissue damage during application
- Reliable, simple and inexpensive manual suction unit for clearing airway obstructions
- Simple disposal eliminates cross contamination concerns



Rev 08262020